


# SUPRAMID

Black , pseudo-monofilament or monofilament , non absorbable polyamide surgical suture  
 Polyamide chir. Nahtmaterial. Schwarz pseudo-monofil oder monofil, nicht resorbierbar



**Reverse Cutting**  
 Schneidende Nadel (DS)

**3/8 CIRCLE**  
 KREIS 3/8



**Reverse Cutting**  
 Schneidende Nadel (HS)

**1/2 CIRCLE**  
 KREIS 1/2

## ▼ 12mm



METRIC	USP	SUTURE LENGTH FADEN-LÄNGE	CODE ARTIKEL
0.7	6/0	75cm	<b>S620-4</b>
1.0	5/0	75cm	<b>S611</b>

## ▼ 18mm



METRIC	USP	SUTURE LENGTH FADEN-LÄNGE	CODE ARTIKEL
2.0	3/0	75cm	<b>S523</b>

## ▼ 16mm



METRIC	USP	SUTURE LENGTH FADEN-LÄNGE	CODE ARTIKEL
1.0	5/0	75cm	<b>S621</b>
1.5	4/0	75cm	<b>S622</b>
2.0	3/0	75cm	<b>S623</b>

## ▼ 19mm



METRIC	USP	SUTURE LENGTH FADEN-LÄNGE	CODE ARTIKEL
1.0	5/0	75cm	<b>S631</b>
1.5	4/0	75cm	<b>S632</b>
2.0	3/0	75cm	<b>S633</b>

\* available upon request / auf Anfrage erhältlich