


ABSORBABLE  
SUTURES  
NON ABSORBABLE  
SUTURES

# Dental Quirurgics



medipac<sup>®</sup>  
manufacturing

Investing in Innovation & High Quality




Medipac S.A.  
Industrial Area PO. Box 1,  
611 00, Kifissia, Greece  
T: +30 23410 71991  
F: +30 23410 71979  
E: info@medipac.gr  
W: www.medipac.gr

# ABSORBABLE SUTURES

# Dental Quirurgics

# NON ABSORBABLE SUTURES




**P.G.A.**  
Braided, violet, synthetic absorbable polyglycolic acid surgical suture, coated

Tensile Strength:		Absorption:
65%	2nd week	60-90 days
35%	3rd week	




**NEOSORB (PGLA)**  
Braided, violet, synthetic absorbable polyglactin 910 surgical suture, coated

Tensile Strength:		Absorption:
50%	3rd week	56-70 days



**NEOSORB RAPID (RPGLA)**  
Braided, undyed, synthetic rapid absorbable polyglactin 910 surgical suture, coated

Tensile Strength:		Absorption:	
50%	6th day	Total lost in:	
		10-14 days	42 days




**MONOSORB**  
Monofilament, violet, synthetic absorbable poly (propionone) surgical suture-PCC

Tensile Strength:		Absorption:	
65% - 70%	4th week	180-220 days	
50% - 60%	6th week		




**MONOFAST**  
Monofilament, violet, synthetic absorbable polyglycolic acid - co - poly caproactone - surgical suture - PGCL


Tensile Strength:		Absorption:	
68% - 78%	1st week	90 - 110 days	
30%	2nd week		



**PROFIMED**  
Undyed monofilament non absorbable, pttle suture



**POLYAMID**  
Blue monofilament, non absorbable polyamide surgical suture (Nylon)




**SILK**  
Black, braided non absorbable wax coated, silk surgical suture



**SUPRAMID**  
Black, Pseudo-monofilament, non absorbable polyamide (Nylon) sterile surgical suture



**PROPYLEN**  
Blue monofilament, non absorbable polypropylen surgical suture



**POLYESTER**  
Green, braided, non absorbable, wax coated polyester surgical suture