

ABSORBABLE SUTURES

Dental Quirurgics

NON ABSORBABLE SUTURES



P.G.A.
Braided, violet, synthetic absorbable polyglycolic acid suture, coated

| Tensile Strength: | Absorption: |
|-------------------|---------------|
| 65% 2nd week | 60-90 days |
| 35% 3rd week | |



NEOSORB (PGLA)
Braided, violet, synthetic absorbable polyglactin 910 surgical suture, coated

| Tensile Strength: | Absorption: |
|-------------------|---------------|
| 50% 3rd week | 56-70 days |



NEOSORB RAPID (RPGLA)
Braided, undyed, synthetic rapid absorbable polyglactin 910 surgical suture, coated

| Tensile Strength: | Absorption: |
|-------------------|---------------------------------|
| 50% 5th day | Total lost in: 10-14 days |
| | 42 days |



MONOSORB
Monofilament, violet, synthetic absorbable poly (p-dioxanone) surgical suture-PDO


| Tensile Strength: | Absorption: |
|----------------------|-----------------|
| 65% -70% 4th week | 180-220 days |
| 50% -60% 6th week | |



MONOFAST
Monofilament, violet, synthetic absorbable polyglycolic acid - co - polycaprolactone surgical suture - PGCL

| Tensile Strength: | Absorption: |
|----------------------|------------------------------|
| 68% -78% 1st week | Total lost in: 28 days |
| 30% 2nd week | 90 -110 days |

PROFIMED
Undyed monofilament non absorbable, ptfe suture



POLYAMID
Blue monofilament non absorbable polyamide surgical suture (nylon)



SILK
Black, braided non absorbable wax coated, silk surgical suture



SUPRAMID
Black, Pseudo-monofilament, non absorbable polyamide (nylon) sterile surgical suture



PROPYLEN
Blue monofilament, non absorbable polypropylene surgical suture



POLYESTER
Green, braided, non absorbable, wax coated polyester surgical suture

